

Bionix OtoClear ear irrigation tips

Product website: <https://bionix.com/otoclear-ear-irrigation.html>



The



OtoClear Tip is designed with 3 holes at the tip which direct the irrigating fluid towards the ear canal walls, and not directly on to the eardrum to reduce risk of any damage. The tip has 3 exit portals for effluent to be drained from the ear.

Use for flushing ear canals as per CARPA manual

Attaches to luer lock syringe

Single use

Supplier

Territory Surgical Supplies (local Darwin supplier)

Product Code EAR7200

<https://tss-nt.com/>